# Information file Euro-pharmat

## Medical Device

### STERICUP SINGLE®

Updated on: 09/02/2020

1. Administrative information about the company			
1.1	<b>Name:</b> APOTHICOM DISTRIBUTION S.A.S.		
1.2	<table border="0"> <tr> <td><b>Full Address:</b> 52 Avenue Edison 75013 Paris France</td> <td> <b>Tel:</b> +33 (0)1 53 61 18 41  <b>Fax:</b> +33 (0)153 61 04 49  <b>E-mail:</b> <a href="mailto:infos@apothicom.org">infos@apothicom.org</a>  <b>Website:</b> <a href="http://www.apothicom.org">www.apothicom.org</a> </td> </tr> </table>	<b>Full Address:</b> 52 Avenue Edison 75013 Paris France	<b>Tel:</b> +33 (0)1 53 61 18 41 <b>Fax:</b> +33 (0)153 61 04 49 <b>E-mail:</b> <a href="mailto:infos@apothicom.org">infos@apothicom.org</a> <b>Website:</b> <a href="http://www.apothicom.org">www.apothicom.org</a>
<b>Full Address:</b> 52 Avenue Edison 75013 Paris France	<b>Tel:</b> +33 (0)1 53 61 18 41 <b>Fax:</b> +33 (0)153 61 04 49 <b>E-mail:</b> <a href="mailto:infos@apothicom.org">infos@apothicom.org</a> <b>Website:</b> <a href="http://www.apothicom.org">www.apothicom.org</a>		
1.3	<table border="0"> <tr> <td><b>Contact details of the person in charge of materiovigilance:</b> Méliné NEGRONI</td> <td> <b>Tel:</b> +33 (0)1 53 61 18 41  <b>Email :</b> <a href="mailto:meline.negronei@apothicom.org">meline.negronei@apothicom.org</a> </td> </tr> </table>	<b>Contact details of the person in charge of materiovigilance:</b> Méliné NEGRONI	<b>Tel:</b> +33 (0)1 53 61 18 41 <b>Email :</b> <a href="mailto:meline.negronei@apothicom.org">meline.negronei@apothicom.org</a>
<b>Contact details of the person in charge of materiovigilance:</b> Méliné NEGRONI	<b>Tel:</b> +33 (0)1 53 61 18 41 <b>Email :</b> <a href="mailto:meline.negronei@apothicom.org">meline.negronei@apothicom.org</a>		

2. Information about the device / equipment	
2.1	<b>Generic Designation:</b> Preparation cooker
2.2	<b>Trade Name:</b> Stericup single
2.3	<b>GMDN Nomenclature Code:</b> 1370
2.4	<b>LPPR Code* (if applicable):</b> N/A
2.5	<p><b>Medical Device Class :</b> Is</p> <p><b>EU applicable Directive:</b> Directive_93/42/CEE    <b>Following Appendix nbr:</b> V, Point 3.</p> <p><b>Notified Body accreditation number:</b> LNE/ G-MED (0459)</p> <p><b>First introduction into EU Market:</b> August 2019</p> <p><b>MD Manufacturer:</b> SO.F.A.P.                      KM14,GP9,ZONE D'ACTIVITE SIDI DAOUD 2046 LA MARSA, TUNIS - TUNISIE</p>
2.6	<p><b>Device Description:</b></p> <p>Stericup® is a sterile and single use cooker with a capacity of 2.5 ml.</p> 

2.7	<b>Description of catalogue references and packaging :</b>		
	<b>Trade Name</b>	<b>References</b>	<b>Description</b>
	Stericup single	AP01-113-99/R-SO	Colour 2.5 ml preparation cooker Overmoulded handle.
	Stericup single	AP01-113-99/SO	Colour 2.5 ml preparation cooker Overmoulded handle.
			<b>Packaging</b>
			Rigid individual sterile blister package Carton of 1 000 units.
			Medium individual sterile blister package Carton of 1 000 units.
			Soft individual sterile blister package Carton of 1 000 units.

2.8	<b>Composition of the device and accessories:</b>	
	<b>Elements:</b>	<b>Materials:</b>
	Preparation cooker	Aluminium
	Overmoulded Handle	Polypropylene
	<b>Additional information about these elements:</b>	
	Principal: the Stericup : Latex free Phthalates (DHP) free PVC free Substance of animal or biological origin free	
	Package: PVC free for the soft blister package (reference AP01-113-99/S-SO)	
	<b>Associated device and accessories to list:</b> N/A	

2.9	<p><b>Purpose of use:</b> The Stericup is intended for people who inject drugs (hereinafter PWID). It has been developed exclusively within the framework of Harm Reduction. It aims at reducing risks/harm associated with the injection of street drugs and medicine diverted from their intended use: viral risk (hepatitis B and C, HIV) and bacterial risks. The “cooker” is a container designed to prepare, and possibly to heat, the solution to be injected. The Stericup has a capacity of 2.5 ml and enables to heat easily a preparation up to 1.5 ml. The overmoulded handle protects fingers from heat and stabilizes the cooker.</p>
-----	---

<b>3. Sterilisation Process:</b>	
<b>Sterile MD:</b>	YES
<b>Device sterilisation mode:</b>	Ethylene Oxide

<b>4. Storage and preservation conditions</b>	
<b>Expiration period:</b>	5 years
Normal conditions of preservation and storage	

<b>5. Safety of use</b>	
5.1	<b>Technical safety:</b> N/A

5.2	<b>Biological safety</b> : N/A
-----	--------------------------------

6. Instructions for use	
6.1	<b>Instructions for use:</b> IFU (Instructions for Use), appendix 1
6.2	<b>Indications:</b> Idem 2-9
6.3	<b>Precautions for use:</b> IFU (Instructions for Use), appendix 1
6.4	<b>Contraindications:</b> N/A

7. Additional information on the device	
	<b><u>Bibliography, clinical trials reports, pharmaco-economic studies, improvement of the provided service, specific recommendations for use</u></b> (treatment restrictions, technical platform, operator qualification etc.): N/A

8. Appendices list (if any)	
	Appendix: IFU

EN

**Stericup®, Sterimix®, Maxicup® and Maximix® are sterile single use cookers.** Their aim is:  
 > To encourage, by its colour, personal identification strategies and reduce the risk of viral transmission (HIV, hepatitis B and C), due to the reuse and sharing of cookers, filters and post-injection pads.  
 > To promote personal and single use of sterile cookers in order to reduce the risks associated with the preparation of injections, namely the risks of bacterial contamination (abscesses, endocarditis...) associated with the use of cookers, filters and post-injection pads.

**STERICUP® OR MAXICUP®?**  
 Stericup® is intended to prepare solutions of up to 1.5 mL. Maxicup® can be used for drugs that need a volume up to 3.5 mL to dissolve. To reduce the risk of group preparation, it is recommended to use Stericup® when its volume is sufficient for one preparation.

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
EN

## INSTRUCTIONS FOR USE

All tools are for personal and single use

1. Open the blister and grab the Stericup® or Maxicup®, making sure not to touch the inside of the cup.

2. **IF YOU USE A STERIFILT®**  
 If you use a syringe without attached needle, the Sterifilt® goes where the needle goes.  
 On all syringes, press and turn slightly the Sterifilt® to snap it into place and ensure a perfect seal.  
 Hold the syringe and Sterifilt® at a slight angle to avoid a suction effect. The plunger is to be pulled back slowly.

3. **IF YOU DON'T USE A STERIFILT®**  
 If the blister includes a cotton filter, use it. Place the needle on the cotton filter and suck up the solution.

4. After the injection, compress the injection site with the provided post-injection pad to encourage healing and prevent infections.

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