Chromatic™ Detection

Chromogenic medium for the enumeration and identification of microorganisms directly from clinical and nonclinical specimens.

DESCRIPTION

Chromatic™ Detection is a chromogenic medium used for the enumeration and identification of microorganisms directly from clinical and nonclinical specimens.

The medium allows to carry out indole test for confirmation of *Escherichia coli*.

TYPICAL FORMULA (g/l)

Peptone	14.0
Tryptone	6.0
Yeast Extract	3.0
Sodium Chloride	5.0
Chromogenic Mix	13.1
Agar	15.0

Final pH 7.2 ± 0.2 at 25°C

METHOD PRINCIPLE

Peptone and tryptone provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Sodium chloride maintains the osmotic balance of the medium. Chromogenic mix allows to identify microorganisms on the basis of the color and morphology of the colonies. Agar is the solidifying agent.

PREPARATION

Dehydrated medium Suspend 56.1 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.

Medium in bottles Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

TEST PROCEDURE

Inoculate the medium by direct streaking or spread plating. Incubate aerobically at 35 ± 2°C for 18-24 hours.

INTERPRETING RESULTS

After incubation observe the color and the morphology of the colonies and interpret the results as indicated in the ID table.

ID Table.

Microorganism	Typical colony color
<i>E. coli</i>	Pink-reddish-mauve
<i>Klebsiella</i> spp, <i>Enterobacter</i> spp, <i>Serratia</i> spp	Green-blue
<i>Proteus</i> spp	Brown
<i>Pseudomonas</i> spp	Yellowish-green
<i>S. aureus</i>	Cream
<i>E. faecalis</i>	Green-turquoise

See pictures in Appendix I.

APPEARANCE

Dehydrated medium: fine, dry, homogeneous, free of extraneous material, beige.

Prepared medium: slightly opalescent, amber.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 2 years.
Medium in bottles: 1 year.
Ready-to-use plates: 4 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.
Inoculum for productivity: 50-100 CFU
Incubation conditions: aerobically at 35 ± 2°C for 18-24 hours.

QC Table.

Microorganism		Growth	Specification
<i>Escherichia coli</i>	ATCC® 25922	Good	Pink-reddish-mauve colonies
<i>Klebsiella pneumoniae</i>	ATCC® 13883	Good	Green-blue colonies
<i>Proteus mirabilis</i>	ATCC® 12453	Good	Brown colonies
<i>Pseudomonas aeruginosa</i>	ATCC® 27853	Good	Yellowish-green colonies
<i>Staphylococcus aureus</i>	ATCC® 25923	Good	Cream colonies
<i>Enterococcus faecalis</i>	ATCC® 29212	Good	Green-turquoise colonies

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

- Merlino J., S. Siarakas, G.J. Robertson, G.R. Funnel, T. Gottlieb, and R. Bradbury (1996) Evaluation of Colorex Orientation for differentiation and presumptive identification of Gram-negative bacilli and Enterococcus species. J. Clin. Microbiol. 34: 1788-1793.
- Samra Z., et al. (1998) Evaluation of use of a new chromogenic Agar in detection of urinary tract pathogens. J. Clin. Microbiol. 36: 990-994.

PRESENTATION

		Contents	Ref.
Chromatic™ Detection	90 mm ready-to-use plates	20 plates	11611
Chromatic™ Detection	Bottles	6 x 100 ml bottles	481130
Chromatic™ Detection	Dehydrated medium	500 g of powder	610612
Chromatic™ Detection	Dehydrated medium	100 g of powder	620612
Chromatic™ Detection	Dehydrated medium	5 kg of powder	6106125

TABLE OF SYMBOLS

LOT Batch code	 Keep away from sunlight	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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