
	<b>TITLE: DEVICE MASTER RECORD</b>	
	<b>QUALITY FILE</b>	<b>Ref. no.: QF 04-5</b>
	<b>Revision: 1</b>	<b>Page 1 of 1</b>
	<b>Issue date: 27.09.2022</b>	<b>Effective date: 27.09.2022</b>

**SUCTION SET**

REF. NO. **1128 3527 80 ML**



**COMPONENT DESCRIPTION**

**COMPONENT NO.**

**1. SUCTION PIECE**

rigid tube ID/OD 3.0/4.0, length: 110 mm, bent  
material: PVC

0.18011040

**2. HANDLE SHORT**

handle white ID 4.0 mm, 160 mm length with **vacuum control**  
material: Styrolux (SBC)

4.21016040

**3. TUBE**

Multi-Lumen transparent tube ID/OD 5.6/8.0 mm,  
material: PVC 73°ShoreA, DEHP – FREE,  
• 1128 3527 80 ML – length 3400 mm

0.56340080

**4. PROXIMAL END**

funnel (funnel) 8.0 PVC  
material: PVC 50°ShoreA, DEHP – FREE, blue

0.00000080

**PACKAGING MATERIAL**

Innerbag 180 x 380 mm

2.18038004

Peelpack (packing machine) 200 x 350 mm

5.20035410

Outerbag 310 x 600 mm

**25 pcs/bag**

2.31060010

Box 386 x 386 x 260 mm

**50 pcs/box**

6.38638626

Label on box 85 x 131 mm

7.00085131

Work instruction

QI 01

Packaging instruction

QI 04

Comments:

Glued

Not glued

Suction piece and handle

X

**Suction piece and handle short put into inner bag.**