





Product data sheet - 10.20.64

26.01.2016

General Data			
Included in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), 450.00.45 (Neuro), 450.00.47 (Cardio), 450.00.53 (Cardio neu), 450.00.54 (OMF Instr), 450.00.58 (Cranialis), 450.00.59 (Spinalis), <b>Single sheets</b> 451.06.85 (Brustdiss),		
Last update	2008-11-06		
Title	<b>Needle Holder</b>		
Article number	<b>10.20.64</b>	Author	Baby-Crile-Wood "HM"
Description	Needle Holder, delicate pattern, Tungsten Carbide inserts, gold-plated handles		
Length	15 cm - 6"	Weight	
		Diameter	
Properties	straight, fine cross serrated, F 5-0, 6-0	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826016098		
GMDN	12726		
Risk category EC	1		
UMDNS	12-726		
Net weight	0.0290		



III CMF III SURGICAL III NEURO+SPINE


**medicon**<sup>®</sup>

Product data sheet - 10.20.64

26.01.2016

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---