





Product data sheet - 49.40.35

24.11.2015

General Data			
Included in	Catalogues 450.00.25 (Kat. 1), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.40 (Kat. 5), Brochures 451.02.30 (Koshima), 451.02.60 (Container), 451.02.61 (Container), 451.03.30 (Servotronic), 451.03.60 (Osteotome), 451.08.50 (Mikro Säg), Single sheets 451.06.06 (Koshima), 451.06.40 (Zylinder-O),		
Parts list	49.46.35 (1), 49.47.01 (1),		
Last update	2014-08-25		
Title	Aluminium Sterilization Container, complete		
Article number	49.40.35	Author	
Description	Aluminium Sterilization Container, complete: Botton silver, Colour of lid: silver, incl. paper filter. Extremely effective drying phase of the container interior, silicone sealings ensure long-term sterility, coloured lids and identification labels for better organization.		
Length		Weight	
		Diameter	
Properties	Length a: 300 mm, Width b: 135 mm, Height c: 47 mm, Lid perforated - Bottom non-perforated	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	C
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPU01		
Operation manuals	G494035_2000_02		
GTIN	4046826051334		
GMDN	13730		
Risk category EC	1		



III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 49.40.35

24.11.2015

UMDNS	13-730
Net weight	0.6400



III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 49.40.35

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	3.0255_3.3315		
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material	1.4301	X5CrNi18-10	Cr 18,0 - 19,5% Ni 8,0 - 10,5%
Others parts material			
Coating material	ELOXIERT	Künstlich erzeugte Oxidschicht bei Aluminium mit und ohne Farbzusatz	

Delivery mode	
---------------	--