





Product data sheet - 16.15.18

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu),		
Last update	2011-05-03		
Title	Hemostatic Forceps		
Article number	16.15.18	Author	Olivecrona (Coller)
Description	Hemostatic Forceps		
Length	18 cm - 7 1/8	Weight	
		Diameter	
Properties	curved, delicate pattern	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826024420		
GMDN	34948		
Risk category EC	1		
UMDNS	10-865		
Net weight	0.0440		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 16.15.18

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---