





Product data sheet - 52.81.00

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu),		
Last update	2008-01-24		
Title	Uterine Curette		
Article number	52.81.00	Author	Recamier
Description	Uterine Curette		
Length	31 cm - 12 1/8"	Weight	
		Diameter	
Properties	sharp, rigid	Tip / Working end	5 mm wide
Article link		Sales unit	
Accessories		Number / Size	fig. 00
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826069087		
GMDN	32594		
Risk category EC	1		
UMDNS	11-084		
Net weight	0.0620		



III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 52.81.00

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material	1.4301	X5CrNi18-10	Cr 18,0 - 19,5% Ni 8,0 - 10,5%
Others parts material			
Coating material			

Delivery mode

