





Product data sheet - 66.02.02

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.48 (Plastic), 450.00.57 (Plastic), 450.00.60 (Plastic neu),		
Last update	2008-01-23		
Title	Nasal Speculum		
Article number	66.02.02	Author	Hartmann
Description	Nasal Speculum		
Length	15 cm - 6	Weight	
		Diameter	
Properties		Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	# 2
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826109820		
GMDN	35350		
Risk category EC	1		
UMDNS	13-664		
Net weight	0.0780		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 66.02.02

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---