





Product data sheet - 202.21.11

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), 450.00.47 (Cardio), 450.00.48 (Plastic), 450.00.53 (Cardio neu), 450.00.57 (Plastic), 450.00.60 (Plastic neu), <b>Brochures</b> 451.03.23 (Medicut),		
Last update	2008-01-24		
Title	<b>Iris Scissors</b>		
Article number	<b>202.21.11</b>	Author	supercut
Description	Iris Scissors		
Length	11.5 cm - 4 1/2	Weight	
		Diameter	
Properties	curved	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826171193		
GMDN	13488		
Risk category EC	1		
UMDNS	13-488		
Net weight	0.0200		




III CMF III SURGICAL III NEURO+SPINE

**medicon**<sup>®</sup>

Product data sheet - 202.21.11

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4117	X38CrMoV15	C 0,35 - 0,40% Cr 14 - 15%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material	T_GOLD	Teil des Instrumentes ist Vergoldet Goldschicht galvanisch max 0,004mm dick	

Delivery mode	
---------------	--