





Product data sheet - 06.05.18

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5),		
Last update	2008-01-24		
Title	Tissue Forceps		
Article number	06.05.18	Author	Standard
Description	Tissue Forceps		
Length	18 cm - 7 1/8	Weight	
		Diameter	
Properties	straight	Tip / Working end	1 x 2 teeth
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826006488		
GMDN	11785		
Risk category EC	1		
UMDNS	14-257		
Net weight	0.0420		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 06.05.18

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4024	X15Cr13	C 0,12 - 0,17% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---