





Product data sheet - 06.02.14

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), 450.00.46 (Sinus), 450.00.54 (OMF Instr),		
Last update	2008-01-24		
Title	<b>Dressing Forceps</b>		
Article number	<b>06.02.14</b>	Author	
Description	Dressing Forceps		
Length	14.5 cm - 5 3/4	Weight	
		Diameter	
Properties	straight, slim	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826006259		
GMDN	11777		
Risk category EC	1		
UMDNS	14-257		
Net weight	0.0250		



III CMF III SURGICAL III NEURO+SPINE


**medicon**<sup>®</sup>

Product data sheet - 06.02.14

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4024	X15Cr13	C 0,12 - 0,17% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---