





Product data sheet - 16.11.18

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.47 (Cardio), 450.00.53 (Cardio neu),		
Last update	2011-10-27		
Title	<b>Hemostatic Forceps</b>		
Article number	<b>16.11.18</b>	Author	Rochester-Pean
Description	Hemostatic Forceps		
Length	18 cm	Weight	
		Diameter	
Properties	curved	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826024017		
GMDN	34948		
Risk category EC	1		
UMDNS	10-865		
Net weight	0.0500		



III CMF III SURGICAL III NEURO+SPINE


**medicon**<sup>®</sup>

Product data sheet - 16.11.18

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---