





Product data sheet - 02.20.10

24.11.2015

| General Data | | | |
|---|---|---|----------|
| Included in | Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), 450.00.48 (Plastic), 450.00.57 (Plastic), 450.00.60 (Plastic neu), | | |
| Last update | 2012-06-13 | | |
| Title | Iris Scissors | | |
| Article number | 02.20.10 | Author | standard |
| Description | Iris Scissors | | |
| Length | 10.5 cm - 4 1/8 | Weight | |
| | | Diameter | |
| Properties | straight | Tip / Working end | |
| Article link | | Sales unit | |
| Accessories | | Number / Size | |
| Complete picture | | Detail images | |
|  | |  | |
| Quality management data | | | |
| Preparation instructions | MPK01 | | |
| Operation manuals | <i>No operation manuals available.</i> | | |
| GTIN | 4046826001513 | | |
| GMDN | 13488 | | |
| Risk category EC | 1 | | |
| UMDNS | 13-488 | | |
| Net weight | 0.0140 | | |




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 02.20.10

24.11.2015

| Materials | | | |
|-------------------------|-----------------|-------------|--------------------------------|
| | Material number | Description | Elements |
| Main component material | 1.4117 | X38CrMoV15 | C 0,35 - 0,40% Cr 14 - 15% |
| Turned part material | 1.4305 | X8CrNiS18-9 | Cr 17,0 - 19,0% Ni 8,0 - 10,0% |
| Handle material | | | |
| Others parts material | | | |
| Coating material | | | |

| | |
|---------------|---|
| Delivery mode |  |
|---------------|---|