





Product data sheet - 29.51.33

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), 450.00.50 (Ortho), 450.00.54 (OMF Instr), 450.00.63 (Ortho),		
Last update	2011-11-23		
Title	Sharp Curette, oval		
Article number	29.51.33	Author	Volkmann
Description	Sharp Curette, oval		
Length	17 cm - 6 3/4	Weight	
		Diameter	
Properties		Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	# 000
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826041342		
GMDN	31335		
Risk category EC	1		
UMDNS	11-084		
Net weight	0.0490		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 29.51.33

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material	1.4301	X5CrNi18-10	Cr 18,0 - 19,5% Ni 8,0 - 10,5%
Others parts material			
Coating material			

Delivery mode	
---------------	---