




Product data sheet - 52.05.12

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu),		
Last update	2008-01-24		
Title	<b>Vaginal Speculum</b>		
Article number	<b>52.05.12</b>	Author	Pederson
Description	Vaginal Speculum		
Length		Weight	
		Diameter	
Properties	large	Tip / Working end	120 x 25 mm
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826063900		
GMDN	35352		
Risk category EC	1		
UMDNS	13-666		
Net weight	0.1870		



III CMF III SURGICAL III NEURO+SPINE

**medicon**<sup>®</sup>

Product data sheet - 52.05.12

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4301	X5CrNi18-10	Cr 18,0 - 19,5% Ni 8,0 - 10,5%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material			

Delivery mode

