





Product data sheet - 52.44.50

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu),		
Last update	2011-04-19		
Title	<b>Tenaculum Forceps</b>		
Article number	<b>52.44.50</b>	Author	Pozzi
Description	Tenaculum Forceps		
Length	25 cm - 10	Weight	
		Diameter	
Properties	straight, screw lock	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826067274		
GMDN	13998		
Risk category EC	1		
UMDNS	16-453		
Net weight	0.0970		



III CMF III SURGICAL III NEURO+SPINE

**medicon**<sup>®</sup>

Product data sheet - 52.44.50

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material			

Delivery mode

