




Product data sheet - 68.24.35

24.11.2015

General Data			
Included in	Catalogues 450.00.30 (Kat. 8), 450.00.48 (Plastic), 450.00.54 (OMF Instr), 450.00.57 (Plastic), 450.00.60 (Plastic neu), Brochures 451.02.36 (Cayman), 451.03.30 (Servotronic), 451.03.31 (Servo II), 451.08.13 (Distractions), 451.08.50 (Mikro Säg), 451.09.01 (bimax), 451.09.02 (bimax), 451.09.10 (Mainzer), 451.09.31 (Verona), 451.09.50 (Maastricht),		
Last update	2015-11-11		
Title	Saw blade		
Article number	68.24.35	Author	
Description	Saw blade for Micro Osseoscalpel Saw according to Sachse		
Length		Weight	
		Diameter	
Properties	Hight of saw blade: 12 mm, thickness of blade: 0.35mm	Tip / Working end	
Article link	68.24.05, 68.24.06	Sales unit	Packaging unit 12 pieces
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826119256		
GMDN	45518		
Risk category EC	2a		
UMDNS	13-447		
Net weight	0.0070		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 68.24.35

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4310	X10CrNi18-8	C 0,05 - 0,15% Cr 16,0 - 19,0% Ni 6,0 - 9,5%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	-------------------------------------------------------------------------------------