





Product data sheet - 68.84.00

24.11.2015

General Data			
Included in	Catalogues 450.00.48 (Plastic), 450.00.52 (implantArt), 450.00.54 (OMF Instr), 450.00.55 (OMF Impl), 450.00.57 (Plastic), 450.00.60 (Plastic neu), 450.00.61 (implantArt), Brochures 451.09.01 (bimax), 451.09.02 (bimax),		
Last update	2011-03-18		
Title	Reference pin		
Article number	68.84.00	Author	Kretschmer
Description	Reference pin for Le Fort I Osteotomy		
Length	20 mm	Weight	
		Diameter	
Properties	Easy pin insertion by means of a special holder, compatible with all current surgical hand pieces (2,35mm shaft), precise determination of vertical dimensions by a fixed reference point, no irritation of the reference point during surgery (swelling of the skin, etc.), easy and precise final checking of the desired vertical position before and during fixation	Tip / Working end	
Article link		Sales unit	Packing unit 10 pieces
Accessories	68.84.05, 68.84.06	Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	G688400_2002_04		
GTIN	4046826124786		



III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 68.84.00

24.11.2015

GMDN	46640
Risk category EC	2a
UMDNS	12-443
Net weight	0.0100




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 68.84.00

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---