



Product data sheet - 47.20.25

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5),		
Last update	2015-03-03		
Title	Kidney Bowl		
Article number	47.20.25	Author	
Description	Kidney Bowl		
Length		Weight	
		Diameter	
Properties	250 mm, 18/8 stainless steel	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826050672		
GMDN	42893		
Risk category EC	1		
UMDNS	13-692		
Net weight	0.2580		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 47.20.25

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4301	X5CrNi18-10	Cr 18,0 - 19,5% Ni 8,0 - 10,5%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---