



## Product Technical Form

Product Name:	Vacutore® Air/Water Bottle Tubing with CO2 Extend
Compatibility	For Olympus® EVIS 140/160/180/190/240/260/290 Series Endoscopes
Product Code:	GAR025CO2EXT
Manufacturer and Address:	GA Health Company Limited Unit 18, 21/F, Metropole Square, 2 On Yiu Street, Shatin, N.T. HONG KONG
MDD Classification:	Class IIa sterile
Quality System:	ISO 13485
UMDNS Code:	14202
Intended Use:	Air/Water Bottle Tubing is intended to connect an air/CO <sub>2</sub> source, a sterile water source (water bottle), and an endoscope to supply air/CO <sub>2</sub> and water during gastrointestinal endoscopic procedures.
Product Photo:	A clear plastic medical tubing assembly. It features a long, thin clear tube that loops back. One end has a clear plastic connector with a blue stopcock. The other end has a larger clear plastic connector with a blue stopcock. The tubing is shown against a white background.
Material:	Free of latex, animal origin or biological products, and DEHP
Sterilization and standard	EO Sterilization, ISO 11135
Packaging:	Primary Packaging: Paper-film pouch bag, 1 piece per bag Secondary Packaging: Color box, 10 pieces per box Carton box, 4 boxes (40 pieces) per carton



Storage:	Keep away from sunlight and keep dry
Disposal:	After use, this product may be a potential biohazard. Handle and dispose of the product in accordance with accepted medical practice and applicable local, state and federal laws and regulations.