

Product Data Sheet

Orthopaedic Suction Set

Pennine 
One Product. One Patient.

Surgical Suction

Reference	AEEFS1111
Description	Sterile, single use, easyflow filtered suction system, for orthopaedic surgery.
Intended Use	Orthopaedic Suction Set is used as a conduit to remove bodily fluids and secretions, surgical tissue debris and irrigation fluids. They are used in surgical sites and body cavities to facilitate observation and/or to clear an obstruction. (GMDN 35917)
Duration of Use	Normally intended for continuous use of less than 60 minutes, but will be patient specific depending on the duration of the surgical procedure.
Materials	Filter/Chamber/Handle - PC Tubing and connectors - PVC Suction Probes - PVC All materials used are DEHP and LATEX Free
Product Information	Supplied with: <ul style="list-style-type: none"> • 7mm/3m Flexi-Rib Tubing; • 3 Suction Probes - 2 curved, 1 straight; • Spare Filter <p>Suitable for removal of bone fragments during orthopaedic surgery;</p> <p>Eliminates blockages;</p> <p>Rapid filter exchange;</p> <p>Ergonomic handle design ensures comfort and control;</p> <p>Reversible suction tip provides extra length for intra-femoral suction.</p>
Packaging	10 units in the outer carton.
Quality System	ISO 13485
EU Classification	Class IIa
Sterilised by	Ethylene Oxide
Shelf Life	5 years



Tubing Length: 3m
Tubing Diameter: 7mm (ID)
Tubing Shore Hardness: 78A



15cm Curved Probe;
23cm Curved Probe;
23 cm Straight Probe.
Probe Diameter: 7.8mm (OD)

Storage and Safety:

Toxicity	Manufactured from biocompatible materials
Storage	Store in a clean, cool, dry place away from direct sunlight



GCMEDICA ENTERPRISE LTD., (WUXI)
Loujin Industrial Park, Shuofang, Wuxi 214143, Jiangsu, China



Luxus Lebenswelt GmbH
Kochstr.1, 47877, Willich, Germany

© 2020 Pennine Healthcare or one of its subsidiaries. All rights reserved.
The information given in this sheet is not intended to be exhaustive.
The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.



Distributed By: Ivor Shaw Ltd trading as Pennine Healthcare
300 City Gate, London Road, Derby, DE24 8WY, United Kingdom
T: +44 (0) 1332 794880 | www.penninehealthcare.co.uk