

BD SmartSite™ needle-free connector

SKU/REF

2000E7D



SKU/REF Name

BD SmartSite™ needle-free connector

Specifications

### Specifications Table

Length (cm)	2.5
Priming Volume (mL)	0.1
Maximum Pressure Rating	325 PSI@10 mL per second
Gravity Flow Rate	>6L/hr
Needle-Free Connector	BD SmartSite™
Change Interval	7days or 200 activations
Type of Luer	Male Spin Lock
DEHP	DEHP-Free
Sterilization Method	Ethylene Oxide
Country of Manufacture	Mexico
Classification	IIa
Quantity per Box	50
Quantity per Case	100
Shelf Life	3 Years